


CP 98/1024

PA 204903



4

# THE UNITED STATES OF AMERICA

TO ALL TO WHOM THESE PRESENTS SHALL COME;

UNITED STATES DEPARTMENT OF COMMERCE

United States Patent and Trademark Office

February 10, 2000

THIS IS TO CERTIFY THAT ANNEXED HERETO IS A TRUE COPY FROM THE RECORDS OF THE UNITED STATES PATENT AND TRADEMARK OFFICE OF THOSE PAPERS OF THE BELOW IDENTIFIED PATENT APPLICATION THAT MET THE REQUIREMENTS TO BE GRANTED A FILING DATE UNDER 35 USC 111.

APPLICATION NUMBER: 60/107,006

FILING DATE: November 04, 1998

REC'D 29 FEB 2000

WIPO PCT

PRIORITY DOCUMENT  
SUBMITTED OR TRANSMITTED IN  
COMPLIANCE WITH  
RULE 17.1(a) OR (b)



By Authority of the  
COMMISSIONER OF PATENTS AND TRADEMARKS

P. SWAIN  
Certifying Officer

# PROVISIONAL APPLICATION COVER SHEET

This is a request for filing a PROVISIONAL APPLICATION under 37 C.F.R. § 1.53(c)				Docket Number: 018845.0158	Type a plus sign (+) inside this box: +
INVENTOR(S)/APPLICANT(S)					
Last Name	First Name	Middle Initial	Residence (City and Either State or Foreign Country)		
BL	David	N.	1089 Goodson Crescent, Oakville, Ontario L6H 4A7 Canada		
SKEA	Danna	L.	5940 Glen Erin Drive, Mississauga, Ontario L5M 5W9 Canada		
HEDGE	Phyllis	R.	6806 Wellington Road 34, R.R. #22, Cambridge, Ontario N3C 2V4 Canada		
TITLE OF THE INVENTION (280 characters max)					
Methods for the Production of TcRγδ <sup>+</sup> T Cells					
CORRESPONDENCE ADDRESS					
James Remenick Baker & Botts, L.L.P. The Warner, Suite 1300 1299 Pennsylvania Ave., N.W.					
State	Washington, D.C.	Zip Code	20004-2400	Country	USA
ENCLOSED APPLICATION PARTS (check all that apply)					
<input checked="" type="checkbox"/> Specification	Number of Pages	4	<input type="checkbox"/> Small Entity Statement		
<input checked="" type="checkbox"/> Drawings	Number of Sheets	4	<input type="checkbox"/> Other (Specify) _____		
METHOD OF PAYMENT (check one)					
<input type="checkbox"/> A check or money order is enclosed to cover the Provisional filing fees.			PROVISIONAL FILING FEE AMOUNT		
<input type="checkbox"/> The Commissioner is hereby authorized to charge any additional filing fees and credit Deposit Account Number: 02-0375.			\$		

The invention was made by an agency of the United States Government or under a contract with an agency of the United States Government.

☒ No.

☐ Yes, the name of the US Government agency and the Government contract number are:

Respectfully submitted,

SIGNATURE \_\_\_\_\_

DATE: November 4, 1998

TYPED OR PRINTED NAME James Remenick \_\_\_\_\_

REGISTRATION NUMBER (if appropriate): 36,902 \_\_\_\_\_

\_\_\_\_\_  
Additional inventors are being named on separately numbered sheets attached hereto.

PROVISIONAL APPLICATION FILING ONLY

## METHODS FOR THE PRODUCTION OF TcR $\gamma\delta^+$ T CELLS

The following methods relate to the large scale, *ex vivo* expansion of TcR $\gamma\delta^+$  T cells in liquid culture in the absence of a stromal layer. The starting material consists of low density mononuclear cells (LDMNC) from human peripheral blood. The LDMNC may be further fractionated by (1) enrichment for CD4 $^+$  T cells, (2) enrichment for T cells together with depletion of TcR $\alpha\beta^+$  T cells, or (3) not further fractionated. The cells are cultured in medium containing conditioned medium (XLCM<sup>TM</sup>), human sera or plasma (P), concanavalin A (con A), interleukin-2 (IL-2), and/or interleukin-4 (IL-4). At frequent intervals the cells are counted and reseeded with fresh medium, XLCM<sup>TM</sup>, P, con A, IL-2, and/or IL-4. The percent of cells expressing a particular surface marker is determined using specific antibodies and flow cytometry.

### Example 1. CD4e - XLCM<sup>TM</sup>/P --- TcR $\gamma\delta^+$ T cells

Low density mononuclear cells were isolated from adult peripheral blood (LDMNC).

CD4 $^+$  cells were enriched from the LDMNC by negative selection using lineage specific antibodies and immunomagnetic affinity chromatography (CD4c).

The CD4e cells were expanded in culture medium containing 5% XLCM<sup>TM</sup> + 5% P.

The cells expanded more than 100,000-fold in four weeks.

After 21 days, more than 50% of the cultured cells were TcR $\gamma\delta^+$ .

The majority of the TcR $\gamma\delta^+$  T cells were V $\delta^+$ .

See Figure 1.

**Example 2. TeAbd - XLCM™/P ---- TcRγδ<sup>+</sup> T cells**

Low density mononuclear cells were isolated from adult peripheral blood (LDMNC).

T cells were enriched and TcRαβ<sup>+</sup> cells were depleted from LDMNC by negative selection using lineage specific antibodies and immunomagnetic affinity chromatography (TeAbd).

The TeAbd were expanded in culture medium containing 5% XLCM™ + 5% P.

The cells expanded more than 100,000-fold in four weeks.

After 8 days, more than 50% of the cultured cells were TcRγδ<sup>+</sup>.

After 15 days, more than 80% of the cultured cells were TcRγδ<sup>+</sup>.

The majority of the TcRγδ<sup>+</sup> T cells were Vδ<sup>+</sup>.

See Figure 2.

Advantages: The TcRγδ<sup>+</sup> T cells expand more rapidly, expand to greater levels, and are more pure.

**Example 3. LDMNC - XLCM™/P - IL-2/IL-4/P ---- TcRγδ<sup>+</sup> T cells**

Low density mononuclear cells were isolated from adult peripheral blood (LDMNC).

The LDMNC were not further fractionated.

The LDMNC were expanded for 5 days in culture medium containing 5% XLCM™ + 5% P, following which they were divided, and half were continuously cultured in XLCM™/P, while the other half were washed and sub-cultured in 10 ng/ml IL-2 + 10 ng/ml IL-4 + 5% P.

In both cases, the cells expanded more than 100,000-fold in four weeks.

However, the different conditions gave rise to different kinds of cells.

Less than 5% of the cells cultured continuously in XLCM™/P were TcRγδ<sup>+</sup>, while more than 50% of the cells cultured in XLCM™/P then sub-cultured in IL-2/IL-4/P were TcRγδ<sup>+</sup>.

See Figure 3.

Advantages: The starting cell number can be very low since no initial fractionation is required. The sub-culture gets rid of XLCM™ components, e.g., con A, mezerein, other known or unknown factors.

**Example 4. LDMNC → Con A/IL-2/IL-4/P → IL-2/IL-4/P → TcRγδ<sup>+</sup> T cells**

Low density mononuclear cells were isolated from adult peripheral blood (LDMNC).

The LDMNC were not further fractionated.

The LDMNC were expanded for 5 days in culture medium containing 20 ug/ml concanavalin A + 10 ng/ml IL-2 + 10 ng/ml IL-4 + 5% P, following which they were divided, and half were continuously cultured in concanavalin A + 10 ng/ml IL-2 + 10 ng/ml IL-4 + 5% P.

In both cases, the cells expanded more than 100,000-fold in four weeks.

However, the different conditions gave rise to different kinds of cells.

Less than 5% of the cells cultured continuously in XLCM™/P were TcRγδ<sup>+</sup>, while more than 50% of the cells cultured in concanavalin A/IL-2/IL-4/P then sub-cultured in IL-2/IL-4/P were TcRγδ<sup>+</sup>.

See Figure 4.


Advantages: The starting cell number can be very low since no initial fractionation is required. The culture conditions are completely defined (i.e., no conditioned medium).

We Claim

1. A method for the enrichment and culture of  $\gamma\delta^+$  T cells and cells produced by the method.
2. A method for the enrichment and culture of V $\delta$ -1 T cells and/or V $\delta$ -2 T cells and cells produced by the method.

20407006 140400

United States Patent & Trademark Office  
Office of Initial Patent Examination -- Scanning Division



Application deficiencies were found during scanning:

☐ Page(s) \_\_\_\_\_ of No Declaration were not present  
for scanning. (Document title)

☐ Page(s) \_\_\_\_\_ of \_\_\_\_\_ were not present  
for scanning. (Document title)

☐ Scanned copy is best available.

60407006 131103

Figure 1






Figure 2




Figure 3

LDMNC

XLCM™/P followed by IL-2/IL-4/P




Figure 4

LDMNC

Con A/IL-2/IL-4/P followed by IL-2/IL-4/P

